



12 Channel touch screen ECG machine

12channel ECG machine with 8 inch Touch and color screen,Multi-language report formats
Meditech EKG 1212T has large capacity (1000 exam),12 Leads ECG simultaneous acquisition,digital signal processor,Built in battery and AC/DC power supply

Features

8" TFT Color Screen 210 mm Thermal Paper printing Support Chinese, English, Spanish and Turkish Operating interface, can print out Multi-language report formats

Include multiform printing speeds such as 5mm/s, 6.25mm/s ,10mm/s, 12.5mm/s, 25mm/s, 50mm/s. Provide graphic result format for visual diagnostic A4 paper direct print (connect with printer)

It has build-in storage with large capacity,1000 exam 12 Leads ECG simultaneous acquisition, digital signal processor, AC filter, baseline filter and EMG filter of the ECG signals Simultaneous display of 3/6/12 leads ECG waveform, and the state of printing mode, sensitivity, speed, filter etc, to be easy for interpretation Include multiform printing modes and formats, such as 12,6+1(rhythm),6 rhythm 12 rhythm 10 rhythm 8 rhythm 6 Manual Freeze etc; The length of wave printing can be adjusted and the function of periodic printing is included to satisfy different applications The instrument has the function of regular auto-analysis and auto-interpretation of ECG waveform parameter, provide measurement parameters such as HR, P-R interval, P Duration, QRS Duration, T Duration, Q-T interval, Q-Tc, P Axis, QRS Axis, T Axis, R(V5), S(V1), R(V5)+S(V1), etc auto-interpretation conclusion to help reducing doctor's abundance and improve working efficiency

AC/DC power supply includes both

Specifications

