



## BabySono Fetal Doppler home health care

BabySono Pocket Fetal Doppler is a portable and reliable fetal doppler that allows a mother-to-be listen to their unborn child's first heartbeats.

Available in Orange and purple

### ▼ Features

The BabySono comes with an audio output which can be connected with headphone for simple auscultation or simple listening, or be connected to a recorder device to record baby's heartbeat to a CD, MP3, a PC, a cell phone or another storage medium of your choice to share the happiness with your friends and family members.

### Features

Listen to your unborn baby's heartbeat from the comfort of your own home  
Can be used from 12 weeks into the pregnancy  
Record the unborn baby's heartbeat using recording cable  
Record the Mother's heartbeat and play it back to help comfort the new born baby  
Conforms to EC Directive 93/42/EEC(MDD) concerning medical devices  
This product uses a 9V alkaline battery for its power source

### Package Includes

BabySono Fetal Doppler Unit  
Recording Cable  
Cleaning Material  
warranty card

### Option cable

Audio cable for all kind of tablets and Smart phones and workable with most of recorders & Doppler's Apps

### ▼ Specifications

Extended Specifications Ultrasound frequency 2.5MHz Battery Voltage: 9V, IEC6F22 9V alkaline Ultrasonic emitting power 10mW/cm<sup>2</sup> Overall sensitivity at the distances 200mm from the face of the transducer greater than 90dB Spatial-peak temporal-peak pressure less than 8.6kPa Output power less than 26mW Effective area of the ultrasonic transducer active element 5.3cm<sup>2</sup>±20% Audio Output Power less than 0.2W Audio out socket 3.5mm in diameter 16-months

Chargeable battery and Charger:

External Soft light Rechargeable Speaker

Audio cable for all kind of tablets and Smart phones and workable with most of recorders & Doppler's Apps