

Automated external defibrillator AED Defi5cs

Selectable Energy AED , Easy to use, Reliable technology&analysis algorithm, Clear Voice prompt, Clear Voice prompt, Wide-range environments

Features

Programmable Energy Selectable Energy for Adult use or Child Storage Large memory capacity up to 1500 case Self Test Daily,Weekly,Monthly Self test Data Transfer Infrared data transmission Memory Keep track of your data

Specifications

Type: Automatic External Defibrillator Self-test: battery, processor, memory, ECG acquisition system and defibrillator. Mode: Adult, Child Waveform: Biphasic truncated exponential Energy: 200-200 Joules max. Energy Sequence Selectable: 30, 30, 50 J 50, 50, 75 J 150, 150, 200J Charge time: Less than 6 sec. to 150J Less than 8 sec. to 200J Voice Prompt:Extended Visual Indicators: LED prompts Control:Two buttons: ON/OFF, Shock ECG storage: 1,500 events.

Data transmission Infrared